

The logo for APMKorea, featuring the text "APMKorea" in white, bold, sans-serif font centered within a blue rectangular box with a subtle drop shadow.

NIBP OEM Module

(Model : M_NIBP)

Contents

- 1. Concept**
- 2. Composition**
- 3. PCB Dimension**
- 4. Specification**
- 5. Protocol**
- 6. Cuff & Electrical Connection**

Rev V1.00

1. Concept

Oscillometry : The oscillometric method of blood pressure measurement is a non-invasive method that monitors the amplitude of cuff pressure changes during cuff deflation to determine arterial blood pressure. The cuff pressure is first elevated above the patient systolic blood pressure level and the cuff begins to deflate at a certain rate. The initial rise in amplitude of these pressure fluctuations during cuff deflation corresponds closely to the systolic blood pressure. As the cuff is further deflated, these pressure fluctuations increase in amplitude until a peak is reached which is usually referred to as the mean arterial pressure(MAP). As cuff deflation continues, the diastolic pressure can be determined based upon the rapidly diminishing amplitude of the pressure fluctuations. Thus systolic, MAP and diastolic blood pressures can be accurately obtained by supervising the pressure fluctuations while controlling the cuff deflation rate.

2. Composition

2-1. Module Selection Model : M_NIBP

2-2. Components

| Discription | Accessory | | Fig. No. | Order No. |
|----------------------|-----------------------------|---|----------|----------------|
| M_NIBP Components | NIBP Board | | Fig 1 | M_NIBP_PC B |
| | RS232C & Power Harness | Harness to Connect RS-232C Communication & Supply Power | Fig 2 | M_NIBP_C1 |
| | UART TTL & Power Harness | Harness to Connect UART TTL Communication & Supply Power | Fig 3 | M_NIBP_C2 |
| | Cuff | Adult Set(A:Large, B:Medium, C:Small) | Fig 4 | Cuff_Adt_* |
| | | Neonate Set(A:Large, B:Medium, C:Small) | Fig 5 | Cuff_Neo_* |
| | Hose | Silicon Tube 1meter | Fig 6 | D3_1000 |
| | PC Program CD | WindowsXP(only for evaluation test module) | Fig 7 | M_NIBP_SW |

2-3. Evaluation Kit

| Discription | Accessory | | Fig. No. | Order No. | |
|-----------------------------|---------------------------|--|------------|-----------|-------|
| M_NIBP Evaluation Kit | NIBP Board | | Fig 1 | M_NIBP_EV | |
| | RS232C & Power Harness | Harness to Connect RS-232C Communication & Supply Power | Fig 2 | | |
| | Cuff | Adult | Cuff_Adt_B | | Fig 4 |
| | | Neonate | Cuff_Neo_B | | Fig 5 |
| | Hose | Silicon Tube 1meter | | | Fig 6 |
| | PC Program CD | WindowsXP(only for evaluation test module) | | | Fig 7 |

2-4. Figures




Fig 1 NIBP Board(M_NIBP_PCB)

J4 is an Amp part number 103635-9(90° connector).




Fig 1 NIBP Board(M_NIBP_PCB)
J4 is an Molex part number 5268-10(90° connector).




Fig 2 RS232C & Power Harness(M_NIBP_C1)




Fig 3 UART TTL & Power Harness(M_NIBP_C2)




Fig 2 RS232C & Power Harness(M_NIBP_C1)



Fig 4 Adult Cuff(Cuff_Adt_*)



Fig 5 Neonate Cuff(Cuff_neo_*)





Fig 6 Silicon Hose

3. PCB Dimension



| | | | |
|-------------------------------|--|---|--------------|
| Dimension | 105.4 × 69.4 × 38mm(Width ×Length ×Height) | | |
| Weight | Appro. 150g | | |
| Power Source | DC from 6.0 to 15Volt ± 5% | | |
| Power consumption | Maximum12Watt, Typical 6Watt | | |
| Measuring Method | Oscillometric | | |
| Measuring range | SYSTOLIC of Blood Pressure | Adult | 40 – 270mmHg |
| | | Pediatric | 40 – 150mmHg |
| | | Neonate | 40 – 130mmHg |
| | MAP of Blood Pressure | Adult | 30 – 260mmHg |
| | | Pediatric | 30 – 140mmHg |
| | | Neonate | 30 – 120mmHg |
| | DIASTOLIC of Blood Pressure | Adult | 20 – 250mmHg |
| | | Pediatric | 20 – 130mmHg |
| | | Neonate | 20 – 110mmHg |
| | Pulse Rate | Measuring range | 30 - 300bpm |
| Measuring Accuracy | Blood Pressure | ±3mmHg between 0 mmHg and 300 mmHg for operating conditions between 0°C and 50°C. | |
| | Pulse Rate | ± 2% or ± 3 BPM, whichever is greater | |
| Initial Inflation Pressure | Adult | 180mmHg(Default), Variable 120-280mmHg | |
| | Pediatric | 130mmHg(Default), Variable 100-160mmHg | |
| | Neonate | 120mmHg(Default), Variable 80-140mmHg | |
| Startup Initialization Period | | 5 seconds | |



4. Specification

| | | |
|---|---|---|
| Patient Safety | Maximum cuff inflation time is limited to 35 seconds | |
| | cuff pressure exceeds 300 mmHg (Adult & Pediatric modes) or 150mmHg (Neonate mode) at any time | |
| | the cuff has been inflated for 180 seconds (Adult & Pediatric modes) or 90 seconds (Neonate mode) | |
| Operation Condition | 0°C to 50°C, 15% to 95% non-condensing humidity | |
| Storage Condition | -20°C to 65°C, 15% to 95% non-condensing humidity | |
| Altitude | Measurement accuracy is not affected by altitude | |
| Interfacing | UART TTL | Asynchronous Operation Baud rate 9600bps, 8 Data bit, 1 Stop bit, no parity bit |
| | UART RS-232C | |
| The Module meets all relevant parts of the following Safety Standards: <ul style="list-style-type: none"> • IEC60601-1:1997 • IEC/EN60601-2-30:1999/2000 • AAMI SP10:1992/2002 • EN1060-1:1996 • EN1060-3:1997 | | |

5. Protocol

5.1 Protocol

- UART [TTL Level (0~5[V])] / [RS-232C Level]
- Asynchronous Operation
- Baud rate : 9600bps
- 8 Data bit
- 1 Stop bit
- no parity bit



5.2 Host Command Packet Structure(Host ⇒ Module)

| HOST START BYTE | COMMAND BYTE | DATA BYTE(S) | CHECKSUM BYTE |
|-------------------------------|--------------|--|--|
| “.” character (0x3A) 1Byte | 1Byte | Supplemental command data. The number of data bytes varies according to the command. | LSB 1Byte [0x100 - BYTE(Startbyte + Command bytes + Data bytes)] |

COMMAND DATA DESCRIPTION

| COMMAND | DATA | DESCRIPTION |
|---------|-----------------------|---|
| 0x17 | Pressure Data(2Bytes) | Sets initial pressure for next BP Data |
| 0x20 | None | Starts Adult Mode BP |
| 0x87 | None | Starts Pediatric Mode BP |
| 0x28 | None | Starts Neonate Mode BP |
| 0x79 | 0x01, 0x00 | Aborts current BP reading |
| 0x79 | 0x05, 0x00 | Returns current cuff pressure |
| 0x79 | 0x03, 0x00 | Returns data from last BP measurement |
| 0x0C | B0, B1, B2 | External control of pump & valves This command should never be used when the cuff is attached to a patient. B0 = Pump (0x00 = Off, 0x01 = On) B1 = Control Valve (0x00 = Open, 0x01 = Closed) B2 = Dump Valve (0x00 = Open, 0x01 = Closed) |

Sending the Starts BP command will automatically place the module in the Adult monitoring mode, causing excessive cuff pressures on a neonate.

It is the responsibility of the Host system to always determine Adult, Pediatric or Neonate mode when initiating a BP measurement.

5.3 Module Packet Structure(Module ⇒ Host)

| MODULE START BYTE | PACKET BYTE | DATA BYTE(S) | CHECKSUM BYTE |
|----------------------|--------------------------------------|---|--|
| “>” character (0x3E) | total number of bytes in the packet. | the Module data response to the COMMAND issued by the Host. | LSB 1Byte [0x100 - BYTE(Startbyte + Command bytes + Data bytes)] |

5.4 Example to set initial inflate to 180mmHg pressure(Decimal 180 ⇒ Hex 0xB4)

Host ⇒ Module

| START | COMMAND | DATA | CHECKSUM |
|-------|---------|------------|----------|
| 0x3A | 0x17 | 0xB4, 0x00 | 0xFB |

Module ⇒ Host("OK" packets)

| START | PACKET | DATA | CHECKSUM |
|-------|--------|------|----------|
| 0x3E | 0x04 | 0x4F | 0x6F |
| 0x3E | 0x04 | 0x4B | 0x73 |

5.5 Example to start adult mode BP

Host ⇒ Module

| START | COMMAND | DATA | CHECKSUM |
|-------|---------|--------|----------|
| 0x3A | 0x20 | (None) | 0xA6 |

Module ⇒ Host (acknowledges "O" about the start command)

| START | PACKET | DATA | CHECKSUM |
|-------|--------|------|----------|
| 0x3E | 0x04 | 0x4F | 0x6F |

Module ⇒ Host ("K" after BP measurement is completed)

| START | PACKET | DATA | CHECKSUM |
|-------|--------|------|----------|
| 0x3E | 0x04 | 0x4B | 0x73 |

When the Module is taking a BP measurement, the only commands that are "valid" are the abort BP and the returns current cuff pressure commands. The Module will respond to all other commands with a "Busy" data packet (data byte = "B") as shown below.

Module ⇒ Host ("B" module is busy)

| START | PACKET | DATA | CHECKSUM |
|-------|--------|------|----------|
| 0x3E | 0x04 | 0x42 | 0x7C |

5.6 Example to start Pediatric mode BP

Host ⇒ Module

| START | COMMAND | DATA | CHECKSUM |
|-------|---------|--------|----------|
| 0x3A | 0x87 | (None) | 0x3F |

Module ⇒ Host (acknowledges "O" about the start command)

| START | PACKET | DATA | CHECKSUM |
|-------|--------|------|----------|
| 0x3E | 0x04 | 0x4F | 0x6F |

Module ⇒ Host ("K" after BP measurement is completed)

| START | PACKET | DATA | CHECKSUM |
|-------|--------|------|----------|
| 0x3E | 0x04 | 0x4B | 0x73 |

When the Module is taking a BP measurement, the only commands that are "valid" are the abort BP and the returns current cuff pressure commands. The Module will respond to all other commands with a "Busy" data packet (data byte = "B") as shown below.

Module ⇒ Host ("B" module is busy)

| START | PACKET | DATA | CHECKSUM |
|-------|--------|------|----------|
| 0x3E | 0x04 | 0x42 | 0x7C |

5.7 Example to start Neonate mode BP

Host ⇒ Module

| START | COMMAND | DATA | CHECKSUM |
|-------|---------|--------|----------|
| 0x3A | 0x28 | (None) | 0x9E |

Module ⇒ Host (acknowledges "O" about the start command)

| START | PACKET | DATA | CHECKSUM |
|-------|--------|------|----------|
| 0x3E | 0x04 | 0x4F | 0x6F |

Module ⇒ Host ("K" after BP measurement is completed)

| START | PACKET | DATA | CHECKSUM |
|-------|--------|------|----------|
| 0x3E | 0x04 | 0x4B | 0x73 |

When the Module is taking a BP measurement, the only commands that are "valid" are the abort BP and the returns current cuff pressure commands. The Module will respond to all other commands with a "Busy" data packet (data byte = "B") as shown below.

Module ⇒ Host ("B" module is busy)

| START | PACKET | DATA | CHECKSUM |
|-------|--------|------|----------|
| 0x3E | 0x04 | 0x42 | 0x7C |

5.8 Example to abort BP

Host ⇒ Module

| START | COMMAND | DATA | CHECKSUM |
|-------|---------|------------|----------|
| 0x3A | 0x79 | 0x01, 0x00 | 0x4C |

Module ⇒ Host("AK" packets : Module Response when BP measurement in progress)

| START | PACKET | DATA | CHECKSUM |
|-------|--------|------|----------|
| 0x3E | 0x04 | 0x41 | 0x7D |
| 0x3E | 0x04 | 0x4B | 0x73 |

Module ⇒ Host("A" packets : Module Response when BP measurement not in progress)

| START | PACKET | DATA | CHECKSUM |
|-------|--------|------|----------|
| 0x3E | 0x04 | 0x41 | 0x7D |

5.9 Example to return current cuff pressure

Host ⇒ Module

| START | COMMAND | DATA | CHECKSUM |
|-------|---------|------------|----------|
| 0x3A | 0x79 | 0x05, 0x00 | 0x48 |

Module ⇒ Host(current pressure of 258mmHg)

| START | PACKET | DATA | CHECKSUM |
|-------|--------|------------|----------|
| 0x3E | 0x05 | 0x02, 0x01 | 0xBA |

Module ⇒ Host(current pressure of 142mmHg)

| START | PACKET | DATA | CHECKSUM |
|-------|--------|------------|----------|
| 0x3E | 0x05 | 0x8E, 0x00 | 0x2F |

5.10 Example to return data from last BP measurement

Host ⇒ Module

| START | COMMAND | DATA | CHECKSUM |
|-------|---------|------------|----------|
| 0x3A | 0x79 | 0x03, 0x00 | 0x4A |

Module ⇒ Host

| START | PACKET | DATA | CHECKSUM |
|-------|--------|--|--|
| 0x3E | 0x18 | SYS, DIA, N1, N2, N3, HR, MAP, EC, N4, N5 | 0x100 – Modulo 256 sum of all bytes |

DATA PACKET

| Data | Description | Number of Bytes |
|------|--|-----------------|
| SYS | Systolic value in mmHg (unsigned integer, LSB first) | 2 |
| DIA | DDD Diastolic value in mmHg (unsigned integer, LSB first) | 2 |
| N1 | Not used | 1 |
| N2 | Not used | 1 |
| N3 | Not used | 8 |
| HR | Heart Rate in beats per minute (unsigned integer, LSB first) | 2 |
| MAP | Mean Arterial Pressure (MAP) in mmHg (unsigned integer, LSB first) | 2 |
| EC | Error Code (unsigned byte) | 1 |
| N4 | Not used | 1 |
| N5 | Not used | 1 |

Error Code List & Definitions

| Error Code | Description | Corrective Action: |
|------------|---|---|
| 0 = 0x00 | Good BP Reading | None |
| 1 = 0x01 | Weak or no oscillometric signal | Check that the cuff is in the correct position. Check the patient. Check that the cuff is properly tightened. Check that there is no excessive clothing between the arm and the cuff. Check that the correct size cuff is being applied. |
| 2 = 0x02 | Artifact / erratic oscillometric signal | The patient may have been moving too much. Check that the cuff is in the correct position. Check that the correct size cuff is being applied. |
| 4 = 0x04 | Exceeded measurement time limit | The patient may have been moving too much. Check that the cuff is properly tightened. Check that the cuff is in the correct position. Check that the correct size cuff is being applied. Check that there is no excessive clothing between the arm and the cuff. |
| 85 = 0x55 | Pneumatic Blockage | Check that the hose has no sharp bends or is pinched. Check that the patient is not lying on the cuff. Check that the cuff is in the correct position. |
| 86=0x56 | BP reading terminated by user | Check the patient. Take another BP reading. |
| 87 = 0x57 | Inflate Timeout, Air Leak or Loose Cuff | Check that the hose is connected to the system and the cuff. Check that the cuff is properly tightened. Check that the cuff is in the correct position. Check that the correct size cuff is being applied. Check that the cuff is not leaking air. Check that the hose connections are not damaged or loose. |

| | | |
|-----------|---|--|
| 89 = 0x59 | Cuff Overpressure | Check that the correct size cuff is being applied. Check that the hose has no sharp bends or is pinched. Check that the cuff is in the correct position. Check that the patient is not lying on the cuff. |
| 90 = 0x5A | Power supply out of range or other hardware problem | Check that V_{PUMP} and V_{LOGIC} are within the voltage specifications. Check the data/power input connection. Service may be required. Call a manufacturer representative. |
| 97 = 0x61 | Transducer out of range | It is not correct calibration, need to recalibrate. Service may be required. Call a manufacturer representative. |
| 98 = 0x62 | ADC out of range | Service may be required. Call a manufacturer representative. |
| 99 = 0x63 | EEPROM calibration data failure | It is not correct calibration, need to recalibrate. Service may be required. Call a manufacturer representative. |

5.11 CONTROL_PNEUMATICS

Caution : This command should never be used when the cuff is attached to a patient.

Host ⇒ Module

| START | COMMAND | DATA | CHECKSUM |
|-------|---------|------------|-------------------------------------|
| 0x3A | 0x0C | B0, B1, B2 | 0x100 – Modulo 256 sum of all bytes |

Data definitions:

- B0 = Pump (0x00 = Off, 0x01 = On)
- B1 = Control Valve (0x00 = Open, 0x01 = Closed)
- B2 = Dump Valve (0x00 = Open, 0x01 = Closed)

Example to close the Control Valve and the Dump Valve:


Host ⇒ Module

| START | COMMAND | DATA | CHECKSUM |
|-------|---------|------------------|----------|
| 0x3A | 0x0C | 0x00, 0x01, 0x01 | 0xB8 |

Module ⇒ Host("OK" packets)

| START | PACKET | DATA | CHECKSUM |
|-------|--------|------|----------|
| 0x3E | 0x04 | 0x4F | 0x6F |
| 0x3E | 0x04 | 0x4B | 0x73 |

6. Cuff & Electrical Connection



6.1 Power and communication connector (PCB J4)

| Pin NO.(PCB J4) | Description | Name |
|-----------------|--|----------------------|
| 1 | No Connect or Optional RS-232 TxD (data out) | TXD-RS |
| 2 | No Connect or Optional RS-232 RxD (data in) | RXD-RS |
| 3 | Ground of Pump | GND _{PUMP} |
| 4 | Ground of Pump | GND _{PUMP} |
| 5 | +VCC of Pump | V _{PUMP} |
| 6 | +VCC of Pump | V _{PUMP} |
| 7 | Logic Level RxD (data in) or No Connect | RXD-TTL |
| 8 | Logic Level TxD (data out) or No Connect | TXD-TTL |
| 9 | Ground of Logic | GND _{Logic} |
| 10 | +VCC of Logic | V _{Logic} |

6.2 J4 is an Amp part number 103670-9(straight connector) or 103635-9(90° connector). The mating connector is Amp part number 104257-9. Pins for the mating connector are Amp part number 104480-7.


OR

J4 is an Molex part number 5267-10(straight connector) or 5268-10(90° connector). The mating connector is Molex part number 5264-10. Pins for the mating connector are Molex part number 5263PBT9.

6.3 Caution to communication connection

Do not attempt to connect to both Logic Level and RS232 simultaneously.

If use to connect to Logic Level, have to remove R137 on PCB.



Rear View of PCB