


**APMKorea**

## SpO2 Sensor Module ICOM1 / ICOM2

### Features


- Monitors for both neonates and adults.
- Low power consumption
- Easy interfacing optional RS-232C (ICOM1) or USB (ICOM2) communication.
- Compact size 43x23x5.5 mm.
- Digital signal processing enhances performance under low perfusion, patient motion and other signal interference.




Accessory		ICOM1	ICOM2	Part No.
Finger Sensor	Resuable probe	●	●	ICOM#FA-01
Pin Connector for Sensor	PCB Connector for Sensor Harness	●	●	ICOM#C1
Pin Connector for Comm.	PCB Connector for Comm. Harness	●	●	ICOM#C2
Hole Connector for sensor Wire	Harness Wire	●	●	ICOM#C3
Hole Connector for comm. Wire	Harness Wire	●	●	ICOM#C4
Sensor Harness	PCB Connector for Comm. Harness	●	●	ICOM#C5
Comm. USB Harness	Harness for USB		●	ICOM#2US
Comm. UART Harness	Harness for UART	●		ICOM#1UA
RS-232C Adaptor	Adaptor for PC Com.	★(option)		ICOM#C6
PC Program SW	MS Windows Version	★(option)	●	ICOM#C7



Sensor Harness



UART Harness



**Specification**


Dimension	43 x23 x5.5mm(Width xLength xHeight)		
Weight	Appro. 70g		
Power Source	DC5Volt $\pm$ 5%(Option DC3.3V), 100mA		
Power consumption	ICOM#1 : (Max)183mWatt at DC5V, 121mWatt at DC3.3V		
	ICOM#2 : (Max)230mWatt at DC5V, 152mWatt at DC3.3V		
Measuring range	Pulse	Measuring range	25 - 255bpm
		Accuracy	$\pm$ 1bpm
	SpO2	Measuring range	0 - 100%
		Accuracy	$\pm$ 1%(90-100% )
			$\pm$ 2%(80-89% )
			$\pm$ 3%(60-79% )
Perfusion Pulse Amplitude		Unspecified (Under 60%)	
Perfusion Pulse Amplitude		0.2 - 20%	
Respos Time	$\leq$ 4sec, in case of Pulse Rate 80bpm		
Recovery Time	$\leq$ 4sec, in case of Pulse Rate 80bpm		
Parameter Average	Adaptive Pulse Rate. Ex) 10th, incase of 80bpm		
Wide Patient Range	Neonate - Adult		
Interfacing	UART TTL	Asynchronous Operation Baud rate 38400bps 8 Data bit, 1 Stop bit, no parity bit	
	USB COM Port		

**Protocol**

- UART [TTL Level (0~5[V])]
- Asynchronous Operation
- Baud rate : 38400bps
- 8 Data bit
- 1 Stop bit
- no parity bit


● Communication protocol

- STX = 0xFA
- ETX = 0xFB
- Total 11byte Transmit
- STX + Wave(2byte) + Hr(2byte) + Spo2(2byte) + Bar data(1byte)  
+ Wave gain(1byte) + Status(1byte) + ETX




- Wave, Hr and SpO2 are Decimal to Hex data

Ex1) Wave=3500, Hr=125, SpO2=98, Bar data=5, Wave gain=2, Status=3



Ex2) Wave=1234, Hr=85, SpO2=96, Bar data=10, Wave gain=10, Status=0



## -. Communication data range

Wave=0 ~ 5000 (sensor or finger is open = 0x2500)

Hr=0 ~ 255 (sensor or finger is open = 0x9999)

SpO2=0 ~ 100( sensor or finger is open = 0x9999)

Bar data=0~10


Wave gain=0~10

Status= refer to the Status

## ● Status

Bit	Showing	
Bit0	0	Sensor Open
	1	Sensor In
Bit1	0	Finger Open
	1	Finger In
Bit2	0	Reserve
	1	Reserve
Bit3	0	Reserve
	1	Reserve
Bit4	0	
	1	Moving while measuring
Bit5	0	
	1	No signal for a period
Bit6	0	
	1	Pulse beep
Bit7	0	Reserve
	1	Reserve

SpO2 sensor / Power connection



● Power connection (PCB J1)

Pin NO.(PCB J1)	Name	Description
1	VCC	+5Volt
2	VCC	+5Volt
3	RX(Module part)	0~5Volt
4	TX(Module part)	0~5Volt
5	GND	0Volt
6	GND	0Volt

● Sensor connection (PCB J3)

Pin NO.(PCB J3)	Name	Description	D-sub(9pin) connector's NO.
1	- Photodiode	- Photodiode in Probe	9
2	+ Photodiode	+ Photodiode in Probe	5
3	Sensor	Sensor ID	1
4	GND	GND	6
5	IR Photoled	Input IR drive voltage from module	3
6	RED Photoled	Input RED drive voltage from module	2